Mettler-Toledo AG, Analytical, Sonnenbergstrasse 74, CH-8603 Schwerzenbach, Switzerland

Printed on 100% chlorine-free paper, for the sake of our environment. 11/07
Gedruckt auf 100% chlorfrei hergestelltem Papier. Unserer Umwelt zuliebe.
Imprimé sur papier 100% exempt de chlore. Par souci d'écologie.
Stampato al 100% su carta senza cloro, per la salvaguardia dell'ambiente.
Impreso en papel totalmente exento de cloro, para preservar nuestro entorno.
この説明書は、環境に配慮し無塩素漂白紙に印刷されています。
为保护我们的环境，此纸张为100%无氯纸。
Отпечатано на бумаге, изготовленной без применения хлора. Мы заботимся о
сохранении окружающей среды.

InLab® pH combination polymer electrodes
InLab® pH-Polymer-Einsfabmesskette
InLab® Electrodes combinées de pH avec polymère
InLab® Elettrodi combinati pH con polimero
InLab® Electroodos combinados de pH con polímero
InLab® pH 複合ポリマー電極
InLab® 固态电解质复合pH电极
Комбинированные pH-электроды InLab® с полимерным электролитом

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METTLER TOLEDO

Read this first

The METTLER TOLEDO InLab® pH combination electrodes with XEROLYT® polymer electrolyte are intended for accurate pH measurement in the laboratory or field at ambient pressure. The exact specifications of the pH and temperature range are printed on each individual electrode shaft or can be found under www.mt.com/electrodes. These operating instructions cover the following electrode models:

51 343 100	InLab® Expert	51 343 104	InLab® Expert NTC30
51 343 101	InLab® Expert Pro	51 343 105	InLab® Expert Pt1000
51 343 102	InLab® Expert 2m	51 340 288	InLab® 413 SG / 2M
51 343 103	InLab® Expert DIN	51 340 289	InLab® 413 SG / 10M
51 343 165	InLab® Semi-Micro		

Safety precautions

The electrode is to be used only for the purpose stated above. It is advisable to wear safety glasses when handling cleaning agents. The electrodes contain glass parts.

Getting ready

The electrode is delivered ready for use. Please check visually for broken parts. For electrodes with detachable cable connect electrode head with cable connector by twisting clockwise. For cables and electrodes with MultiPin™ connector both marks (white/gray) must be in line before twisting. Connect other end with the pH meter. If the electrode contains an integral temperature sensor, the cable for the temperature sensor must also be connected. Next remove the wetting cap and rinse with distilled water. Eliminate any air bubbles inside the glass membrane space (internal reference electrolyte) by smoothly shaking the unit in the vertical plane (as for a fever thermometer). The electrode is now ready for calibration. Please follow the detailed calibration procedure, described in the manual of your pH meter. After calibration the electrode is ready for measurement.

Practical hint: After rinsing do not rub the glass membrane, since this can increase response time or damage the electrode.

Storage

Rinse or clean electrode after use and replace the wetting cap, partly filled (1/3) with Friscolyt-B™ (order no: 51 340 053) or 3 mol/L KCl (51 340 049).

Warning: Storage in distilled water will drastically decrease the life time of the electrode.

Trouble-shooting

- check all connections;
- temperature sensor must be connected to the corresponding pH meter input;
- for tips go to www.electrodes.net

Slow response/drift

- clean open apertures with warm tap water;
- after dry storage soak for 24 hours in Friscolyt-B™.

Slope not adjustable

- clean open apertures with warm tap water;
- clean connector.

Contamination by ...

Oily/organic liquids – Rinse with acetone or ethanol (flammable);

Warning: Cleaning with strong acids (eg HCl) will decrease the life time;

Each time the electrode is cleaned, it must then be recalibrated.

Disposal

If no relative official regulations are in force, used or defective InLab® sensors and the packaging can be disposed of as ordinary household waste.

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