

## Instructions for use

# HYCON® Agar Strips

## Commercial Size and Packing:

- 1.44253.0050 | Agar Strips TC | 50 pieces
- 1.44226.0040 | Agar Strips TC-γ | 40 pieces
- 1.44240.0050 | Agar Strips TSM | 50 pieces
- 1.44228.0040 | Agar Strips TCI-γ | 40 pieces
- 1.44255.0050 | Agar Strips PEN | 50 pieces
- 1.44109.0040 | Agar Strips PEN-γ | 40 pieces
- 1.44108.0040 | Agar Strips LAC-γ | 40 pieces
- 1.44242.0050 | Agar Strips YM | 50 pieces
- 1.44243.0050 | Agar Strips SDX | 50 pieces
- 1.44244.0040 | Agar Strips SDX-γ | 40 pieces
- 1.44245.0025 | Agar Strips DG18 | 25 pieces
- 1.44102.0025 | Agar Strips S | 25 pieces
- 1.44099.0025 | Agar Strips C | 25 pieces

## Components

The composition and application of agar strips is outlined in the respective Certificate of Analysis. Certificates are available from the website ([www.merckmillipore.com](http://www.merckmillipore.com)).

## Characteristics

Ready-to-use culture media for assessment of airborne microorganisms with HYCON® Microbial Air Samplers, i.e. RCS® High Flow Touch, RCS® High Flow, RCS® Plus, RCS® Plus Ex, RCS® Standard, and RCS® Isolator. The portfolio includes agar strips for determination of the total colony counts as well as yeasts and molds (also in gamma-irradiated, double-wrapped packaging, respectively), and staphylococci and coliform bacteria. Agar strips for the detection of airborne microorganisms are individually sealed, flexible plastic strips filled with different types of media.

## Application

1. Prior to use the agar strip should be equilibrated to room temperature.
2. Open the wrapper approximately at 1/3 by peeling back the plastic seal at the rounded side of the wrapper. Remove the agar strip with the coated side facing downwards.
3. Insert the agar strip into the opening of the rotor, or the impeller drum according to the directions outlined in the user manual of the respective microbial air sampler.
4. Place the instrument into required position, choose the appropriate sample volume and start the air sampling procedure.
5. When sampling is finished, remove the agar strip and place it back into the original wrapper. Seal the wrapper with an adhesive tape or Cover Slides (Art. No. 1.44111.0100).
6. Label the wrapper e.g. with a waterproof pen for identification.

7. Proceed with incubation according to the following conditions:

Determination of	Incubation Time	Temperature
Total count	1-5 days	30°C-35°C
Yeasts and molds	3-7 days	20°C-25°C
Staphylococci, coliforms	1-3 days*	30°C-35°C

\*During incubation of agar strips type S the typical growth pattern for staphylococci (color change of colonies to yellow due to fermentation of mannitol) is usually achieved after 18-24 hours. Prolonged incubation times may yield to higher colony counts, yet growth might be less specific.

### Important Notes

- Prior to use the unopened strip should be subjected to a visual inspection for dehydration or contamination. Dehydrated or contaminated agar strips should be discarded.
- Practice aseptic technique when handling agar strips.
- The coated surface of the agar strips should face down during incubation in order to avoid the formation of satellites by condensing water.
- Country-specific or internal guidelines may indicate incubation times and temperatures that differ from those outlined above.

### Evaluation

After incubation, the colony-forming units (CFUs) are analyzed by visual inspection directly through the sealed wrapper.

### Storage and shelf-life

- The expiration date can be derived from the Certificate of Analysis, the box label, and the tab of each agar strip.
- Agar strips are stored in original packaging (Styrofoam box).
- Agar strips are stored protected from light. In particular type YM, C and S should be stored protected from light to avoid the inactivation of supplements and dyes.
- Agar strips type S and C are stored cooled at 2°C to 15°C. All other agar strips are stored at 2°C to 25°C.
- Examples of possible storage conditions:
  - Room Temperature (20°C to 25°C). Storage under these conditions does not require equilibration prior to use.
  - Refrigerator (2°C to 8°C).
  - Beverage Cooler (8°C to 15°C).

### Important Notes

- To prevent condensation effects avoid frequent changes between temperature conditions.
- Upon storage agar strips should not be placed near heat sources such as refrigerators with heat-emitting condensators.
- Boxes should be stored with the coated side of the agar strip facing downwards (label is readable on boxes with 50 pieces; label is facing downwards on boxes with 25 pieces).

### Disposal

Living microorganisms grown on agar strips must be inactivated prior to disposal. Agar strips can be autoclaved, burned, or treated with biocides. Where stipulated, the country-specific or internal guidelines must be followed.



We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

EMD Millipore and the M logo are registered trademarks of Merck KGaA, Darmstadt, Germany.

Version: June 2013 Lit. No. TN5340ENUS

### More HYCON® Products for Environmental Monitoring

- HYCON® RCS Air Sampler, i.e. RCS® High Flow Touch, RCS® Isolator, RCS® Plus Ex and RCS® Standard for determining the microbial count in the surrounding air.
- Wide range of accessories such as the RCS® Compressed Gas Adapter for microbial monitoring of compressed gasses.
- HYCON® Airborne Particle Counters i.e. APC SmartTouch, ErgoTouch and ErgoTouch Pro 2.
- HYCON® Validation and Qualification Handbooks.
- HYCON® Contact Slides for surface monitoring.

*Distributed by:*



800-826-8302 [nelsonjameson.com](http://nelsonjameson.com)

EMD Millipore Corporation  
290 Concord Road  
Billerica, MA 01821, U.S.A.  
[mibio@emdmillipore.com](mailto:mibio@emdmillipore.com)  
[www.emdmillipore.com/biomonitoring](http://www.emdmillipore.com/biomonitoring)

Find contact information for your country at:  
[www.emdmillipore.com/offices](http://www.emdmillipore.com/offices)  
For Technical Service, please visit:  
[www.emdmillipore.com/techservice](http://www.emdmillipore.com/techservice)



We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

EMD Millipore and the M logo are registered trademarks of Merck KGaA, Darmstadt, Germany.

Version: June 2013 Lit. No. TN5340ENUS