RCS[®] Compressed Gas Adapter and Nozzle Set

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Merck

Merck Life Science Pvt. Ltd.

50 A, 2nd Phase | Ring Road | Peenya | Bangalore – 560 058

Tel: +91 80 3928 2500

Fax: +91 80 2839 6345

India

Merck KGaA

Frankfurter Strasse 250

64293 Darmstadt

Germany

Merck India is part of Merck KGaA

hycon@merckgroup.com

www.your-contacts.com/biomonitoring

MANUAL HISTORY

The following is a manual history of the
RCS® Compressed Gas Adapter
and Nozzle Set
Revision Date: May 2015

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IMPORTANT NOTE

In the following Merck shall mean Merck KGaA
Darmstadt, Germany.

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Contents

1 Introduction	4
2 Functional Principle	5
3 Operating Instructions	6
3.1 General Notes	6
3.2 Sterilization and Decontamination	6
3.3 Installation and Sampling	6
4 Nozzle Set for RCS® Compressed Gas Adapter	8
4.1 Application Range	8
4.2 Mounting Instructions	9
5 Accessories	10
6 Technical Data	11
7 Limited Warranty	12

Scope of Supply RCS® Compressed Gas Adapter and NozzleSet:

- 1 × RCS® Compressed Gas Adapter
 with reduction nozzle 1 bar;
 (optional) 1 × Carrying Case
 with five reduction nozzles
- 1 × Operation Manual
- 1 × Carrying Case

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1 Introduction

In clean rooms and other controlled environments, compressed air and gases are used for the operation of technical devices and for cleaning of the manufactured products. Microbial monitoring of pressured gas lines is therefore, required for the environmental monitoring program of the clean room and ensures the quality of the products manufactured therein.

RCS® Microbial Air Samplers provide a renowned technology for microbial monitoring of surrounding air. RCS® instruments are based on the principle of centrifugal impaction according to Reuter (Reuter Centrifugal Sampler, RCS®). In general, RCS® Microbial Air Samplers are designed to handle air at atmospheric pressures. The RCS® Compressed Gas Adapter is designed to reduce the air velocity of the compressed gas and therefore, apply the compressed gas pressure-less to the RCS® instrument. Together with the Nozzle Set an inlet pressure range of 0.1 Absolute to 7.0 bar Absolute can be handled.

The RCS® Compressed Gas Adapter is used in combination with the RCS® High Flow Touch, RCS® High Flow, and RCS® Isolator Microbial Air Samplers for investigating the microbiological content of pressurized gases according to ISO8753-7. It is routinely applied to areas where specific microbial count limits must be met due to regulations, legal requirements, and individual standards. According to ISO14698, a sample volume of up to 1 m³ (1000 L) should be used for the microbial monitoring of compressed gas.


The RCS® Compressed Gas Adapter allows for a simple yet effective measurement of the microbial quality of compressed air and gas in the following environments:

- Isolators and laminar flow
- Sterile and aseptic production lines in pharmaceutical industries


2 Functional Principle

The RCS[®] Compressed Gas Adapter consists of a diffusor, a reduction nozzle, and a tube connection that allows for the connection of a compressed gas pipeline to the adapter.

The tube connector of the RCS[®] Compressed Gas Adapter has a thread size of R1/4 and is designed as a universal connection. The tube connector can be replaced easily by a connector of a different size, or a quick coupling connection. It is important that the connector is always securely attached to the reduction nozzle – an important component of the RCS[®] Compressed Gas Adapter.

The gas to be analyzed enters the RCS[®] Compressed Gas Adapter through the tube connection. Thereby, the velocity of the volume flow is restricted by the reduction nozzle. The flow cross-section of the air stream is expanded within the diffusor, and the flow velocity is continuously reduced. Following, the gas is used by the attached RCS[®] Microbial Air Sampler almost pressure-free. An excess gas rate throughout the measurement allows for a constant flow rate during sampling. The excess amount of gas is released through exhaust channels and directed over the protection cap of the RCS[®] Microbial Air Sampler.

3 Operating Instructions

3.1 General Notes

IMPORTANT

The equipment needed for evaluating microbial counts of the compressed gas includes the RCS® Compressed Gas Adapter and a RCS® Microbial Air Sampler type RCS® High Flow Touch, RCS® High Flow, and RCS® Isolator. The RCS® Compressed Gas Adapter is standardly delivered with the 1 bar nozzle installed. The optional Nozzle Set (Art. No.1.44235.0001) expands the application area from 1 bar to 0.1 Absolute to 7.0 bar Absolute.

For the investigation of the total microbial counts in compressed gases, the RCS® Microbial Air Sampler is equipped with a rotor containing an agar strip. For the analysis of dry compressed gases, agar strips type TCI-γ are recommended.

Apply aseptic techniques when handling the instrument and the agar strip to ensure that colony counts are derived from the compressed gas only.

Ensure that all air outlet openings of both the RCS® units and the RCS® Compressed Gas Adapter are not obstructed during operation in order to avoid turbulences, or variations in the flow rate.

3.2 Sterilization and Decontamination

Sterilization and decontamination measures should be established for the RCS® Compressed Gas Adapter and the reduction nozzles and carried out routinely. The components of the adapter are made of aluminium and stainless steel. Cleaning and disinfection can be performed by rinsing or spraying with commercially available cleaning and disinfecting agents. The agent should have a neutral pH value (pH: 5 – 8) to protect the anodized aluminium coating of the diffusor. Please ensure to use cleaning agents compatible with aluminium and stainless steel. Alternatively, the adapter can be gas-sterilized or autoclaved at 121°C. After sterilization, it is recommended to store the adapter in a sterile container or bag.

3.3 Installation and Sampling

The assembly of all components has to take place in a microorganism-free environment. It is recommended to sterilize the RCS® Compressed Gas Adapter before use. Prior to microbial sampling of compressed gas, the RCS® Microbial Air Sampler must be equipped with a rotor containing the agar strip.

The RCS® Compressed Gas Adapter is connected to the pressured gas line. Ensure a tight and close fit between the adapter and pressured gas line. If necessary, disinfect the connecting parts.

For assembly of the RCS® Compressed Gas Adapter with the RCS® Microbial Air Sampler, place the RCS® instrument in a horizontal position on a stable surface. Gently push the diffuser of the adapter over the protection cap of the sampling head of the RCS® Microbial Air Sampler allowing it to be secured in this position. Next, open the pressure valve and make sure that the pressure does not exceed the pressure range of the reduction nozzle. As the compressed air begins to flow through the adapter, it will produce a marked hissing noise that is derived from the expansion cylinder. The gas should escape for approximately 1 minute before the RCS® Microbial Air Sampler is started. After completion of air sampling, close the valve of the compressed line and remove the RCS® Microbial Air Sampler from the adapter.

If required, transfer the RCS® Microbial Air Sampler or the rotor containing the agar strip into a sterile container in order to protect the agar strip from secondary contaminations. Remove the agar strip from the rotor and place it back into the original wrapper. Label the wrapper and proceed with the incubation of the agar strip. If multiple measurements are performed for the same compressed gas outlet, do not close the compressed air valve during replacement of agar strip or rotor.

IMPORTANT NOTES

- The diffuser of the RCS® Compressed Gas Adapter contains a small foot at one side. For installations with the RCS® High Flow Microbial Air Sampler, the foot is facing downwards. For installations with the RCS® High Flow Touch Microbial Air Sampler, the foot is facing upwards.
 - According to ISO 14698, a sample volume of 1 m³ can be used for evaluation of compressed gas. Therefore, 1000 L is selected as sampling volume on the RCS® Instrument.
 - If examinations are performed in non-controlled environments, it is recommended to equip the rotor with the agar strip in a sterile environment and transport it in a sterile container. Alternatively, several rotors can be prepared. The RCS® Air Samplers allow the management of up to 10 rotors that are calibrated individually with the instrument.
 - For more information on the operation of the RCS® Microbial Air Sampler, please refer to the respective user manual.
-

4 Nozzle Set for RCS® Compressed Gas Adapter

4.1 Application Range

The Nozzle Set is used in combination with the RCS® Compressed Gas Adapter (Art. No. 1.44257.0001) for microbial monitoring of compressed gases. The set consists of five reduction nozzles that are designed to extend the application range of the RCS® Compressed Gas Adapter from an air inlet pressure of 1 bar to 0.1 Absolute to 7.0 bar Absolute. In total, there are six nozzles available for the RCS® Compressed Gas Adapter. The 1 bar nozzle is installed upon delivery of the adapter.

For microbial monitoring of compressed gas pipelines, select the appropriate nozzle according to the inlet flow rate of the pressured gas within the gas line. Then, install the nozzle according to the following instruction.

Pressure in bar	0.10 – 0.20 bar Absolute	0.20 – 0.40 bar Absolute	0.40 – 0.80 bar Absolute	0.80 – 1.20 bar Absolute	2.00 – 4.00 bar Absolute	5.00 – 7.00 bar Absolute
Nozzle	0.15 bar	0.30 bar	0.60 bar	1.00 bar	3.00 bar	6.00 bar

Pressure characteristics curves are available upon request.

IMPORTANT NOTES

- Selection of the correct nozzle is mandatory to allow for a flow rate of 100 L/min for the RCS[®] Microbial Air Sampler. A flow rate lower than 100 L/min may force the RCS[®] instrument to pull air from the surrounding environment in addition to the air received from the gas line. This may lead to incorrect results.
- If the pressure of the gas line is at the upper or lower limits of a nozzle, please select the nozzle for the lower pressure range. For example, if the inlet pressure is 0.40 bar use the 0.30 nozzle. This will ensure that the flow rate is maintained and an excess of gas prevents that air is taken in from the environment.


4.2 Mounting Instructions

Required Products:

- 1 × RCS[®] Compressed Gas Adapter
- 1 × Carrying Case with 5 Nozzles
- 1 × User Manual


Required Tools:

- 1 × Jaw wrench WS14 (WS = wrench size)
- 1 × Jaw wrench WS17

Picture 1. Compressed Gas Adapter with replacement nozzle and tools

1. Obtain 2 small adjustable wrenches, or a jaw wrench WS 14 and a jaw wrench WS 17. Place the jaw wrench WS 14 on the nozzle mount and the jaw wrench WS 17 on the barbed hose fitting.
2. Hold the nozzle mount while slightly twisting the jaw wrench on the barbed fitting to remove this piece.
3. Remove the nozzle that is installed inside the cone and replace it with the nozzle chosen with help of the table above. Make sure to install the copper sealing ring on top of the nozzle to ensure a tight fit.
4. Screw the barbed fitting securing the nozzle into place by using the jaw wrenches.

Picture 2. Disassembled tube connector and sealing ring

5 Accessories

Determination of the microbial count in compressed gases requires a RCS® Microbial Air Sampler, RCS® Compressed Gas Adapter, optional the Nozzle Set, as well as, agar strips.

	Article No.	Description
RCS® High Flow Touch	1.44194.0001	Portable air sampler with a flow rate of 100 L/min. Including power supply, serial RS232 cable, USB adapter, RCS® Management Software, rotor, protection cap, carrying case, calibration certificate, quick start guide and user manual
RCS® Isolator	1.44236.0001	Air sampler for fixed installations with a flow rate of 100 L/min. Including sampling unit, operation unit, connection cable with a secured plug, infrared remote control, rotor, protection cap, carrying case and user manual
RCS® Compressed Gas Adapter	1.44257.0001	Autoclave-able adapter for microbial monitoring of compressed gas; designed for a pressure of 1 bar
Nozzle Set for RCS® Compressed Gas Adapter	1.44235.0001	Set of five nozzles to extend the air inlet pressure from 1 bar to 0.1 Absolute – 7.0 bar Absolute
Agar Strips, Type TCI-γ	1.44228.0040	Gamma-irradiated Tryptic Soy Agar with neutralizers, double-wrapped for determination of total count in aseptic environments and in peroxide-containing air. For examinations in aseptic environments. Suited for monitoring of dry gases

6 Technical Data

Principle:	Adapter for the microbial monitoring of compressed gas with the RCS® High Flow Touch, RCS® High Flow and RCS® Isolator Microbial Air Sampler
Range of application:	Inlet pressure 1 bar + 10%. The optional Nozzle Set for the RCS® Compressed Gas Adapter extends the application range from 1 bar to 0.1 Absolute – 7.0 bar Absolute. Higher pressures must be reduced by pressure regulators. The adapter is only suitable for non-flammable and non-toxic gases
Sampling volume:	1000 L (recommended)
Flow rate:	With primary pressure of 1 bar, the flow rate is approx. 180 L/min
Diffusor and Reduction Nozzle:	Anodized aluminum and stainless steel. Autoclave-able and H ₂ O ₂ -sterilizable. Can be disinfected with commercially available detergents (pH: 5 – 8)

7 Limited Warranty

Products of Merck are sold under the terms and conditions of sale of Merck and its affiliates, which may vary from country to country. The products sold hereunder are warranted to be free from defects in workmanship and material at the time of delivery to the customer. Merck makes no warranty or representation, either expressed or implied, with respect to the suitability of a product for a particular purpose. There are no warranties, expressed or implied, which extend beyond the technical specifications of the products.

Merck declares its instruments to be free from defects of materials and workmanship. The warranty period is one (1) year commencing on the date of commissioning of the instrument as evidenced by means of an invoice or other documents.

The warranty covers the rectification of all faults which occur as a proven direct result of faulty materials or workmanship during the warranty period.

The following items are not covered by the warranty:

- 1 Parts subject to normal wear and tear.
- 2 Malfunction caused by accidental or intentional damage, or the use of parts from other manufacturers, or through repairs or modifications carried out by persons not authorized by Merck to carry out such work.
- 3 Defects which occur as a result of normal wear and tear or from using the equipment under conditions not specified in the instructions for use.

The user is not entitled to any remedy in law towards Merck, howsoever arising, apart from those remedies mentioned here. Merck accepts no other liability or warranty towards the user. The legal rights of the purchaser towards the seller remain unaffected by this.

Please be advised that the RCS® instruments are specified for use with HYCON® Agar Media and other HYCON® consumables only. Instrument usage with agar strips or consumables not produced, distributed or authorized by Merck will void any warranty as correct operation of the instrument cannot be guaranteed. Furthermore, Merck will not be liable for any faults or defects originated from using inappropriate agar media or consumables.

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We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose. The M mark is a trademark of Merck KGaA, Darmstadt, Germany.

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