-  (EN) **Molecular Detection Assay 2 - *Listeria***
-  (FR) **Kit de détection moléculaire *Listeria* sp version 2**
-  (DE) **Molekulare Detektion 2 – *Listerien* Nachweis**
-  (IT) **Analisi molecolare per il rilevamento della *Listeria***
-  (ES) **Ensayo de Detección Molecular 2 para *Listeria***
-  (NL) **Moleculair detectie assay 2 *Listeria***
-  (SV) **Molecular Detection Assay 2 – *Listeria***
-  (DA) **Molekylær Detektions Analyse 2 *Listeria***
-  (NO) **Molekylær deteksjonstest 2 for *Listeria***
-  (FI) **Molekyläärinen testisetti 2 - *Listeria***
-  (PT) **Ensaio para Detecção Molecular de *Listeria* 2**
-  (EL) **Δοκιμασία Μοριακής Ανίχνευσης *Listeria* 2**
-  (PL) **Molekularny test do wykrywania 2 – *Listeria***
-  (RU) **Молекулярный анализ 2 - *Listeria***
-  (TR) **Moleküler Tayin Testi *Listeri*-2**
-  (JA) **分子検出アッセイ2 リステリア属菌**
-  (ZH) **李斯特菌分子检测分析 2 - 李斯特菌 (Simplified)**
-  (ZH) **分子檢測套組 2 - 李斯特菌 (Traditional)**
-  (TH) **ชุดทดสอบเชื้อลิสทีเรียระดับโมเลกุล2**
-  (KO) **분자 검출 키트2 - 리스테리아용**
-  (ID) **Deteksi Molekuler Untuk Pengujian *Listeria* 2**

2

Product Instructions

Molecular Detection Assay 2 - *Listeria*

MDA2LIS96

PRODUCT DESCRIPTION AND INTENDED USE

3M™ Molecular Detection Assay 2 - *Listeria* is used with the 3M™ Molecular Detection System for the rapid and specific detection of *Listeria* species in enriched food and environmental samples.

The 3M Molecular Detection Assays use loop-mediated isothermal amplification to rapidly amplify nucleic acid sequences with high specificity and sensitivity, combined with bioluminescence to detect the amplification. Presumptive positive results are reported in real-time while negative results are displayed after the assay is completed. Presumptive positive results should be confirmed using your preferred method or as specified by local regulations^(1, 2, 3).

The 3M Molecular Detection Assay 2 - *Listeria* is intended for use in a laboratory environment by professionals trained in laboratory techniques. 3M has not documented the use of this product in industries other than food or beverage. For example, 3M has not documented this product for testing water, pharmaceutical, cosmetics, clinical or veterinary samples. The 3M Molecular Detection Assay 2 - *Listeria* has not been evaluated with all possible testing protocols or with all possible strains of bacteria.

As with all test methods, the source, formulation and quality of enrichment medium can influence the results. Factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may also influence results. 3M recommends evaluation of the method including enrichment medium, in the user's environment using a sufficient number of samples with particular foods and microbial challenges to ensure that the method meets the user's criteria.

3M has evaluated the 3M Molecular Detection Assay 2 - *Listeria* with Demi-Fraser Broth containing Ferric Ammonium Citrate. A typical formulation of this medium follows below.

Demi-Fraser Broth Base Typical Formula (g/L)

Sodium Chloride	20 g
Sodium Phosphate, dibasic, anhydrous *	9.6 g
Beef Extract	5.0 g
Pancreatic Digest of Casein	5.0 g
Peptic Digest of Animal Tissue	5.0 g
Yeast Extract	5.0 g
Lithium Chloride	3.0 g
Potassium Phosphate, monobasic	1.35 g
Esculin	1.0 g
Acriflavin HCl	0.0125 g
Nalidixic Acid	0.01 g
* Substitute: Sodium Phosphate, dibasic, dihydrate	12.0 g

Fraser Broth Supplement

(Ingredients per 10 mL vial. One vial is added to one liter of basal medium.)

Ferric Ammonium Citrate 0.5 g/10 mL

Final pH 7.2 ±0.2 at 25°C

The 3M™ Molecular Detection Instrument is intended for use with samples that have undergone heat treatment during the assay lysis step, which is designed to destroy organisms present in the sample. Samples that have not been properly heat treated during the assay lysis step may be considered a potential biohazard and should NOT be inserted into the 3M Molecular Detection Instrument.

3M Food Safety is certified to ISO (International Organization for Standardization) 9001 for design and manufacturing.

The 3M Molecular Detection Assay 2 - *Listeria* test kit contains 96 tests, described in Table 1.

Table 1. Kit Components

Item	Identification	Quantity	Contents	Comments
Lysis Solution (LS) tubes	Pink solution in clear tubes	96 (12 strips of 8 tubes)	580 µL of LS per tube	Racked and ready to use
<i>Listeria</i> Reagent tubes	Blue tubes	96 (12 strips of 8 tubes)	Lyophilized specific amplification and detection mix	Ready to use
Extra caps	Blue caps	96 (12 strips of 8 caps)		Ready to use
Reagent Control (RC)	Clear flip-top tubes	16 (2 pouches of 8 individual tubes)	Lyophilized control DNA, amplification and detection mix	Ready to use
Quick Start Guide		1		

The Negative Control, not provided in the kit, is sterile enrichment medium, e.g., Demi-Fraser Broth. Do not use water as a Negative Control.

SAFETY

The user should read, understand and follow all safety information in the instructions for the 3M Molecular Detection System and the 3M Molecular Detection Assay 2 - *Listeria*. Retain the safety instructions for future reference.

⚠ **WARNING:** Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

⚠ **CAUTION:** Indicates a hazardous situation, which, if not avoided, could result in minor or moderate injury and/or property damage.

NOTICE: Indicates a potentially hazardous situation which, if not avoided, could result in property damage.

⚠ WARNING

Do not use the 3M Molecular Detection Assay 2 - *Listeria* in the diagnosis of conditions in humans or animals.

The 3M Molecular Detection Assay 2 - *Listeria* method may generate *Listeria monocytogenes* to levels sufficient to cause stillbirths and fatalities in pregnant women and the immunocompromised, if exposed.

The user must train its personnel in current proper testing techniques: for example, Good Laboratory Practices, ISO 17025⁽⁴⁾, or ISO 7218⁽⁵⁾.

To reduce the risks associated with a false-negative result leading to the release of contaminated product:

- Follow the protocol and perform the tests exactly as stated in the product instructions.
- Store the 3M Molecular Detection Assay 2 - *Listeria* as indicated on the package and in the product instructions.
- Always use the 3M Molecular Detection Assay 2 - *Listeria* by the expiration date.
- Use the 3M Molecular Detection Assay 2 - *Listeria* for food and environmental samples that have been validated internally or by a third party.
- Use the 3M Molecular Detection Assay 2 - *Listeria* only for surfaces, sanitizers, protocols and bacterial strains that have been validated internally or by a third party.
- For an environmental sample containing Neutralizing Buffer with aryl sulfonate complex, perform a 1:2 dilution before testing (1 part sample into 1 part sterile enrichment broth). 3M™ sample handling products which include neutralizing buffer: BPPFV10NB, RS96010NB, RS9604NB, SSL10NB, XLSLSSL10NB, HS10NB and HS119510NB.

To reduce the risks associated with exposure to chemicals and biohazards:

- It is strongly recommended that female laboratory staff be informed of the risk to a developing fetus resulting from infection of the mother through exposure to *Listeria monocytogenes*.
- Perform pathogen testing in a properly equipped laboratory under the control of trained personnel.
- Always follow standard laboratory safety practices, including wearing appropriate protective apparel and eye protection while handling reagents and contaminated samples.
- Avoid contact with the contents of the enrichment media and reagent tubes after amplification.
- Dispose of enriched samples according to current industry standards.

To reduce the risks associated with cross-contamination while preparing the assay:

- Always wear gloves (to protect the user and prevent introduction of nucleases).

To reduce the risks associated with environmental contamination:

- Follow current industry standards for disposal of contaminated waste.

⚠ CAUTION

- Do not exceed the recommended temperature setting on heater.
- Do not exceed the recommended heating time.
- Use an appropriate, calibrated thermometer to verify the 3M™ Molecular Detection Heat Block Insert temperature (e.g., a partial immersion thermometer or digital thermocouple thermometer, not a total immersion thermometer.) The thermometer must be placed in the designated location in the 3M Molecular Detection Heat Block Insert.

NOTICE

To reduce the risks associated with cross-contamination while preparing the assay:

- Use of sterile, aerosol barrier (filtered), molecular biology grade pipette tips is recommended.
- Use a new pipette tip for each sample transfer.
- Use Good Laboratory Practices to transfer the sample from the enrichment to the lysis tube. To avoid pipettor contamination, the user may choose to add an intermediate transfer step. For example, the user can transfer each enriched sample into a sterile tube.
- Use a molecular biology workstation containing germicidal lamp where available.

To reduce the risks associated with a false-positive result:

- Never open tubes post amplification.
- Always dispose of the contaminated tubes by soaking in a 1-5% (v:v in water) household bleach solution for 1 hour and away from the assay preparation area.

Consult the Safety Data Sheet for additional information and local regulations for disposal.

If you have questions about specific applications or procedures, please visit our website at www.3M.com/foodsafety or contact your local 3M representative or distributor.

LIMITATION OF WARRANTIES / LIMITED REMEDY

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, 3M DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any 3M Food Safety Product is defective, 3M or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify 3M within sixty days of discovery of any suspected defects in a product and return it to 3M. Please call Customer Service (1-800-328-1671 in the U.S.) or your official 3M Food Safety representative for a Returned Goods Authorization.

LIMITATION OF 3M LIABILITY

3M WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall 3M's liability under any legal theory exceed the purchase price of the product alleged to be defective.

USER RESPONSIBILITY

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.3M.com/foodsafety, or contact your local 3M representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any 3M Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

To help customers evaluate the method for various food matrices, 3M has developed the 3M™ Molecular Detection Matrix Control kit. When needed, use the Matrix Control (MC) to determine if the matrix has the ability to impact the 3M Molecular Detection Assay 2 - *Listeria* results. Test several samples representative of the matrix, i.e. samples obtained from different origin, during any validation period when adopting the 3M method or when testing new or unknown matrices or matrices that have undergone raw material or process changes.

A matrix can be defined as a type of product with intrinsic properties such as composition and process. Differences between matrices may be as simple as the effects caused by differences in their processing or presentation for example, raw vs. pasteurized; fresh vs. dried, etc.

STORAGE AND DISPOSAL

Store the 3M Molecular Detection Assay 2 - *Listeria* at 2-8°C. Do not freeze. Keep kit away from light during storage. After opening the kit, check that the foil pouch is undamaged. If the pouch is damaged, do not use. After opening, unused reagent tubes should always be stored in the re-sealable pouch with the desiccant inside to maintain stability of the lyophilized reagents. Store resealed pouches at 2-8°C for no longer than 60 days.

Do not use 3M Molecular Detection Assay 2 - *Listeria* past the expiration date. Expiration date and lot number are noted on the outside label of the box. After use, the enrichment medium and the 3M Molecular Detection Assay 2 - *Listeria* tubes can potentially contain pathogenic materials. When testing is complete, follow current industry standards for the disposal of contaminated waste. Consult the Safety Data Sheet for additional information and local regulations for disposal.

INSTRUCTIONS FOR USE

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Periodically decontaminate laboratory benches and equipment (pipettes, cap/decap tools, etc.) with a 1- 5% (v:v in water) household bleach solution or DNA removal solution.

SAMPLE ENRICHMENT

Table 2 presents guidance for the enrichment of food and environmental samples. It is the user's responsibility to validate alternate sampling protocols or dilution ratios to ensure this test method meets the user's criteria.

Foods

1. Allow the Demi-Fraser Broth enrichment medium (includes ferric ammonium citrate) to equilibrate to ambient laboratory temperature.
2. Aseptically combine the enrichment medium and sample according to Table 2. For all meat and highly particulate samples, the use of filter bags is recommended.
3. Homogenize thoroughly by blending, stomaching, or hand mixing for 2 ± 0.2 minutes. Incubate at $37 \pm 1^\circ\text{C}$ according to Table 2.
4. For raw dairy products, transfer 0.1 mL of the primary enrichment into 10 mL of Fraser Broth. Incubate at $37 \pm 1^\circ\text{C}$ for 20-24 hours.

Environmental samples

Sample collection devices can be a sponge hydrated with a neutralizing solution to inactivate the effects of the sanitizers. 3M recommends the use of a biocide-free cellulose sponge. Neutralizing solution can be Dey-Engley (D/E) Neutralizing Broth or Lethen broth. It is recommended to sanitize the area after sampling.

WARNING: Should you select to use Neutralizing Buffer (NB) that contains aryl sulfonate complex as the hydrating solution for the sponge, it is required to perform a 1:2 dilution (1 part sample into 1 part sterile enrichment broth) of the enriched environmental sample before testing in order to reduce the risks associated with a false-negative result leading to the release of contaminated product.

The recommended size of the sampling area to verify the presence or absence of the pathogen on the surface is at least 100 cm^2 (10 cm x 10 cm or 4"x4"). When sampling with a sponge, cover the entire area going in two directions (left to right then up and down) or collect environmental samples following your current sampling protocol or according to the FDA BAM⁽¹⁾, USDA FSIS MLG⁽²⁾ or ISO 18593⁽⁶⁾ guidelines.

1. Allow the Demi-Fraser Broth enrichment medium (includes ferric ammonium citrate) to equilibrate to ambient laboratory temperature.
2. Aseptically combine the enrichment medium and sample, according to Table 2.
3. Homogenize thoroughly by blending, stomaching, or hand mixing for 2 ± 0.2 minutes. Incubate at $37 \pm 1^\circ\text{C}$ for 24-30 hours.

Table 2: Enrichment protocols using Demi-Fraser Broth Enrichment

Sample Matrix	Sample Size	Enrichment Broth Volume (mL)	Enrichment Temperature (°C)	Enrichment Time (hr)				
Heat-processed, cooked, cured meats, poultry, seafood and fish Heat-processed / pasteurized dairy products Produce and vegetables Multi-component foods	25 g	225	37	24-30				
Environmental samples	1 sponge	100 or 225	37	24-30				
	1 swab	10	37	24-30				
Raw meat, poultry, seafood, fish	25 g	475	37	28-32				
Sample Matrix	Primary Enrichment (Demi-Fraser Broth)				Secondary Enrichment (Fraser Broth)			Sample Analysis Volume (a)
	Sample Size	Enrichment Broth Volume (mL)	Enrichment Temperature (°C)	Enrichment Time (hr)	Sample Size	Enrichment Temperature (°C)	Enrichment Time (hr)	
Raw dairy products	25 g	225	37	20-24	Transfer 0.1 mL into 10 mL Fraser Broth	37	20-24	10 µL

(a) Volume of sample transferred to Lysis Solution tubes. Refer to Step 4.6 of Lysis section.

PREPARATION OF THE 3M™ MOLECULAR DETECTION SPEED LOADER TRAY

1. Wet a cloth with a 1-5% (v:v in water) household bleach solution and wipe the 3M™ Molecular Detection Speed Loader Tray.
2. Rinse the 3M Molecular Detection Speed Loader Tray with water.
3. Use a disposable towel to wipe the 3M Molecular Detection Speed Loader Tray dry.
4. Ensure the 3M Molecular Detection Speed Loader Tray is dry before use.

PREPARATION OF THE 3M™ MOLECULAR DETECTION CHILL BLOCK INSERT

Place the 3M™ Molecular Detection Chill Block directly on the laboratory bench; (the 3M™ Molecular Detection Chill Block Tray is not used). Use the chill block at ambient laboratory temperature (20-25°C).

PREPARATION OF THE 3M™ MOLECULAR DETECTION HEAT BLOCK INSERT

Place the 3M™ Molecular Detection Heat Block Insert in a dry double block heater unit. Turn on the dry block heater unit and set the temperature to allow the 3M Molecular Detection Heat Block Insert to reach and maintain a temperature of 100 ±1°C.

NOTE: Depending on the heater unit, allow approximately 30 minutes for the 3M Molecular Detection Heat Block Insert to reach temperature. Using an appropriate, calibrated thermometer (e.g., a partial immersion thermometer or digital thermocouple thermometer, not a total immersion thermometer) placed in the designated location, verify that the 3M Molecular Detection Heat Block Insert is at 100 ±1°C.

PREPARATION OF THE 3M™ MOLECULAR DETECTION INSTRUMENT

1. Launch the 3M™ Molecular Detection Software and log in.
2. Turn on the 3M Molecular Detection Instrument.
3. Create or edit a run with data for each sample. Refer to the 3M Molecular Detection System User Manual for details.


NOTE: The 3M Molecular Detection Instrument must reach and maintain temperature of 60°C before inserting the 3M Molecular Detection Speed Loader Tray with reaction tubes. This heating step takes approximately 20 min and is indicated by an ORANGE light on the instrument's status bar. When the instrument is ready to start a run, the status bar will turn GREEN.

LYSIS


1. Allow the lysis solution (LS) tubes to warm up by setting the rack at room temperature (20-25 °C) overnight (16-18 hours). Alternatives to equilibrate the LS tubes to room temperature are to set the LS tubes on the laboratory bench for at least 2 hours, incubate the LS tubes in a $37 \pm 1^\circ\text{C}$ incubator for 1 hour or place them in a dry double block heater for 30 seconds at 100°C .
2. Invert the capped tubes to mix. Proceed to next step within 4 hrs.
3. Remove the enrichment broth from the incubator.
4. One LS tube is required for each sample and the Negative Control (NC) (sterile enrichment medium) sample.
 - 4.1 LS tube strips can be cut to desired LS tube number. Select the number of individual LS tubes or 8-tube strips needed. Place the LS tubes in an empty rack.
 - 4.2 To avoid cross-contamination, decap one LS tubes strip at a time and use a new pipette tip for each transfer step.
 - 4.3 Transfer enriched sample to LS tubes as described below:

Transfer each enriched sample into an individual LS tube **first**. Transfer the NC **last**.

- 4.4 Use the 3M™ Molecular Detection Cap/Decap Tool-Lysis to decap one LS tube strip - one strip at a time.
- 4.5 Discard the LS tube cap – if lysate will be retained for retest, place the caps into a clean container for re-application after lysis. For processing of retained lysate, see Appendix A.
- 4.6 Transfer 20 μL of sample into a LS tube unless otherwise indicated in Protocol Table 2. e.g. Raw dairy products use 10 μL .
5. Repeat step 4.2 until each individual sample has been added to a corresponding LS tube in the strip.

6. Repeat steps 4.1 to 4.6 as needed, for the number of samples to be tested.
7. When all samples have been transferred, transfer 20 μL of NC (sterile enrichment medium, e.g., Demi-Fraser Broth) into a LS tube. Do not use water as a NC.
8. Verify that the temperature of the 3M Molecular Detection Heat Block Insert is at $100 \pm 1^\circ\text{C}$.
9. Place the uncovered rack of LS tubes in the 3M Molecular Detection Heat Block Insert and heat for 15 ± 1 minutes. During heating, the LS solution will change from pink (cool) to yellow (hot).
Samples that have not been properly heat treated during the assay lysis step may be considered a potential biohazard and should NOT be inserted into the 3M Molecular Detection Instrument.
10. Remove the uncovered rack of LS tubes from the heating block and allow to cool in the 3M Molecular Detection Chill Block Insert at least 5 minutes and a maximum of 10 minutes. The 3M Molecular Detection Chill Block Insert, used at ambient temperature without the Molecular Detection Chill Block Tray should sit directly on the laboratory bench. When cool, the lysis solution will revert to a pink color.
11. Remove the rack of LS tubes from the 3M Molecular Detection Chill Block Insert.


AMPLIFICATION

1. One Reagent tube is required for each sample and the NC.
 - 1.1 Reagent tube strips can be cut to desired tube number. Select the number of individual Reagent tubes or 8-tube strips needed.
 - 1.2 Place Reagent tubes in an empty rack.
 - 1.3 Avoid disturbing the reagent pellets from the bottom of the tubes.
2. Select 1 Reagent Control (RC) tube and place in rack.

3. To avoid cross-contamination, decap one Reagent tube strip at a time and use a new pipette tip for each transfer step.
4. Transfer lysate to Reagent tubes and RC tube as described below:

Transfer each sample lysate into individual Reagent tubes **first** followed by the NC. Hydrate the RC tube **last**.

5. Use the 3M™ Molecular Detection Cap/Decap Tool-Reagent to decap the Reagent tubes—one Reagent tube strip at a time. Discard cap.
 - 5.1 Transfer 20 µL of Sample lysate from the upper ½ of the liquid (avoid precipitate) in the LS tube into corresponding Reagent tube. Dispense at an angle to avoid disturbing the pellets. Mix by gently pipetting up and down 5 times.
 - 5.2 Repeat step 5.1 until individual Sample lysate has been added to a corresponding Reagent tube in the strip.
 - 5.3 Cover the Reagent tubes with the provided extra caps and use the rounded side of the 3M Molecular Detection Cap/Decap Tool-Reagent to apply pressure in a back and forth motion ensuring that the cap is tightly applied.
 - 5.4 Repeat step 5.1 as needed, for the number of samples to be tested.
 - 5.5 When all sample lysates have been transferred, repeat 4.1 to transfer 20 µL of NC lysate into a Reagent tube.
 - 5.6 Transfer **20 µL of NC lysate into a RC tube**. Dispense at an angle to avoid disturbing the pellets. Mix by gently pipetting up and down 5 times.
6. Load capped tubes into a clean and decontaminated 3M Molecular Detection Speed Loader Tray. Close and latch the 3M Molecular Detection Speed Loader Tray lid.

7. Review and confirm the configured run in the 3M Molecular Detection Software.
8. Click the Start button in the software and select instrument for use. The selected instrument's lid automatically opens.
9. Place the 3M Molecular Detection Speed Loader Tray into the 3M Molecular Detection Instrument and close the lid to start the assay. Results are provided within 75 minutes, although positives may be detected sooner.
10. After the assay is complete, remove the 3M Molecular Detection Speed Loader Tray from the 3M Molecular Detection Instrument and dispose of the tubes by soaking in a 1-5% (v:v in water) household bleach solution for 1 hour and away from the assay preparation area.

NOTICE: To minimize the risk of false positives due to cross-contamination, never open reagent tubes containing amplified DNA. This includes Reagent Control, Reagent, and Matrix Control tubes. Always dispose of sealed reagent tubes by soaking in a 1-5% (v:v in water) household bleach solution for 1 hour and away from the assay preparation area.

RESULTS AND INTERPRETATION

An algorithm interprets the light output curve resulting from the detection of the nucleic acid amplification. Results are analyzed automatically by the software and are color-coded based on the result. A Positive or Negative result is determined by analysis of a number of unique curve parameters. Presumptive positive results are reported in real-time while Negative and Inspect results will be displayed after the run is completed.

Presumptive positive samples should be confirmed as per the laboratory standard operating procedures or by following the appropriate reference method confirmation^(1,2,3), beginning with transfer from the primary enrichment to secondary enrichment broth (if applicable), followed by subsequent plating and confirmation of isolates using appropriate biochemical and serological methods.

NOTE: Even a negative sample will not give a zero reading as the system and 3M Molecular Detection Assay 2 - *Listeria* amplification reagents have a “background” relative light unit (RLU).

In the rare event of any unusual light output, the algorithm labels this as “Inspect.” 3M recommends the user to repeat the assay for any Inspect samples. If the result continues to be Inspect, proceed to confirmation test using your preferred method or as specified by local regulations.

If you have questions about specific applications or procedures, please visit our website at www.3M.com/foodsafety or contact your local 3M representative or distributor.

REFERENCES:

1. US Food and Drug Administration Bacteriological Analysis Manual. Chapter 10: Detection and Enumeration of *Listeria monocytogenes* in Foods. Section C-6. April 2011 Version.
2. US Department of Agriculture (USDA) FSIS Microbiology Laboratory Guidebook 8.08. Isolation and Identification of *Listeria monocytogenes* from Red Meat, Poultry and Egg Products, and Environmental Samples. Effective Date: 6 Nov 2012.
3. ISO 11290-1. Microbiology of Food and Animal Feeding Stuffs – Horizontal Method for the Detection and Enumeration of *Listeria monocytogenes*. Amendment 1, 2004-10-15.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
6. ISO 18593. Microbiology of food and animal feeding stuffs – Horizontal methods for sampling techniques from surfaces using contact plates and swabs.

Appendix A. Protocol Interruption: Storage and re-testing of heat-treated lysates

1. To store a heat-treated lysate, re-cap the lysis tube with a clean cap (see “LYSIS”, 4.5)
2. Store at 4 to 8°C for up to 72 hours.
3. Prepare a stored sample for amplification by inverting 2-3 times to mix.
4. Decap the tubes.
5. Place the mixed lysate tubes on 3M Molecular Detection Heat Block Insert and heat at $100 \pm 1^\circ\text{C}$ for 5 ± 1 minutes.
6. Remove the uncovered rack of LS tubes from the heating block and allow to cool in the 3M Molecular Detection Chill Block Insert at least 5 minutes and a maximum of 10 minutes.
7. Continue the protocol at the ‘Amplification’ section detailed above.

EXPLANATION OF PRODUCT LABEL SYMBOLS

Caution or Warning, see product instructions.

Consult product instructions.

The lot in a box represents the lot number.

The hourglass is followed by a month and year which represent the expiration date.

Storage temperature limitations.

3M Food Safety

3M United States

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-800-328-6553

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800-563-2921

3M Europe and MEA

3M Deutschland GmbH
Carl-Schurz - Strasse 1
D41453 Neuss/Germany
+49-2131-14-3000

3M Latin America

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-954-340-8263

3M Asia Pacific

No 1, Yishun Avenue 7
Singapore, 768923
65-64508869

3M Japan

3M Health Care Limited
6-7-29, Kita-Shinagawa
Shinagawa-ku, Tokyo
141-8684 Japan
81-570-011-321

3M Australia

Bldg A, 1 Rivett Road
North Ryde, NSW 2113
Australia
61 1300 363 878

3M Health Care

2510 Conway Ave
St. Paul, MN 55144-1000 USA
www.3M.com/foodsafety

© 2015, 3M. All rights reserved.
3M is a trademark of 3M. Used under license in Canada.
34-8716-8554-0

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RÉFÉRENCES :

1. US Food and Drug Administration Bacteriological Analysis Manual. Chapter 10: Detection and Enumeration of *Listeria monocytogenes* in Foods. Section C-6. April 2011 Version.
2. US Department of Agriculture (USDA) FSIS Microbiology Laboratory Guidebook 8.08. Isolation and Identification of *Listeria monocytogenes* from Red Meat, Poultry and Egg Products, and Environmental Samples. Effective Date: 6 Nov 2012.
3. ISO 11290-1. Microbiology of Food and Animal Feeding Stuffs – Horizontal Method for the Detection and Enumeration of *Listeria monocytogenes*. Amendement 1, 2004-10-15.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
6. ISO 18593. Microbiology of food and animal feeding stuffs – Horizontal methods for sampling techniques from surfaces using contact plates and swabs.

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REFERENTIES:

1. US Food and Drug Administration Bacteriological Analysis Manual. Chapter 10: Detection and Enumeration of *Listeria monocytogenes* in Foods. Section C-6. April 2011 Version.
2. US Department of Agriculture (USDA) FSIS Microbiology Laboratory Guidebook 8.08. Isolation and Identification of *Listeria monocytogenes* from Red Meat, Poultry and Egg Products, and Environmental Samples. Effective Date: 6 Nov 2012.
3. ISO 11290-1. Microbiology of Food and Animal Feeding Stuff – Horizontal Method for the Detection and Enumeration of *Listeria monocytogenes*. Amendment 1, 2004-10-15.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
6. ISO 18593. Microbiology of food and animal feeding stuffs – Horizontal methods for sampling techniques from surfaces using contact plates and swabs.

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REFERENSER:

1. US Food and Drug Administration Bacteriological Analysis Manual. Chapter 10: Detection and Enumeration of *Listeria monocytogenes* in Foods. Section C-6. April 2011 Version.
2. US Department of Agriculture (USDA) FSIS Microbiology Laboratory Guidebook 8.08. Isolation and Identification of *Listeria monocytogenes* from Red Meat, Poultry and Egg Products, and Environmental Samples. Effective date: 6 Nov 2012.
3. ISO 11290-1. Microbiology of Food and Animal Feeding Stuffs – Horizontal Method for the Detection and Enumeration of *Listeria monocytogenes*. Amendment 1, 2004-10-15.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
5. ISO 7218. Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
6. ISO 18593. Microbiology of food and animal feeding stuffs – Horizontal methods for sampling techniques from surfaces using contact plates and swabs.

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1. US Food and Drug Administration Bacteriological Analysis Manual. Chapter 10: Detection and Enumeration of *Listeria monocytogenes* in Foods. Section C-6. April 2011 Version.
2. US Department of Agriculture (USDA) FSIS Microbiology Laboratory Guidebook 8.08. Isolation and Identification of *Listeria monocytogenes* from Red Meat, Poultry and Egg Products, and Environmental Samples. Effective Date: 6 Nov 2012.
3. ISO 11290-1. Microbiology of Food and Animal Feeding Stuffs - Horizontal Method for the Detection and Enumeration of *Listeria monocytogenes*. Amendment 1, 2004-10-15.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
5. ISO 7218. Microbiology of food and animal feeding stuffs - General rules for microbiological examination.
6. ISO 18593. Microbiology of food and animal feeding stuffs - Horizontal methods for sampling techniques from surfaces using contact plates and swabs.

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