

- 
-  **(EN) Fraser Broth Base**
 -  **(FR) Base Bouillon Fraser**
 -  **(DE) Fraser Basis Bouillon**
 -  **(IT) Brodo di base Fraser**
 -  **(ES) Base de Caldo Fraser**
 -  **(NL) Fraser bouillon (basis)**
 -  **(SV) Fraser buljong**
 -  **(DA) Fraser Bouillon Base**
 -  **(NO) Fraser buljong**
 -  **(FI) Fraser Elatusainejauhe**
 -  **(PT) Fraser Caldo Base**
 -  **(EL) Βάση Ζωμού Fraser**
 -  **(PL) Bulion do frezowania**
 -  **(HU) Fraser Tápleves Alap**
 -  **(CS) Fraserův bujón – báze**
 -  **(RO) Bulion bază Fraser**
 -  **(RU) Бульон Фрейзера**
 -  **(TR) Toz Fraser Broth**
 -  **(JA) フレーザー基礎培地**
 -  **(ZH) Fraser 肉汤基础**
 -  **(TH) เฟรเซอร์บรอตเบส**
 -  **(KO) Fraser Broth Base**

Product Instructions

Fraser Broth Base

DESCRIPTION

Medium used for the selective enrichment of *Listeria* species including *L. monocytogenes*.

APPLICATIONS

For Laboratory Use Only. Safety Data Sheet available.

PREPARATION OF MEDIUM

Suspend 55 g. of the powder in 1.0 L of purified water.

Mix thoroughly.

Autoclave at 121°C for 15 minutes.

Cool to room temperature.

Final pH 7.2 ± 0.2 at 25°C.

ADDITIVE

3M Fraser Broth Supplement (5% Ferric Ammonium Citrate, **REF:** BP0220010) is required for differential enrichment of *Listeria*. If the addition of 3M Fraser Broth Supplement is necessary, cool the broth base to room temperature and aseptically add 10 mL of the 3M Fraser Broth Supplement. Mix thoroughly.

STORAGE

2-30°C. Hygroscopic. Keep container tightly closed.

DISPOSAL

Follow established laboratory practices.

Disclaimer: 3M disclaims all express and implied warranties including warranties of merchantability or fitness for a particular use. If product is defective, the exclusive remedy is, at 3M's option, replacement, or refund. Except where prohibited by law, 3M will not be liable for any further loss or damage arising from use of this product.

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.3M.com/foodsafety, or contact your local 3M representative or distributor for more information.

EXPLANATION OF SYMBOLS

www.3M.com/foodsafety/symbols

Distributed by:

NELSON JAMESON INC.

800-826-8302 nelsonjameson.com