

3M Science.
Applied to Life.™

Scrub dot technology to disrupt biofilm
and enhance sample collection

Hydrated with Wide
Spectrum Neutralizer

A new standard of sampling.

Do more with the 3M™ Environmental Scrub Sampler.

Now you can add a sampling solution to your environmental monitoring program that is designed for compatibility with 3M™ Petrifilm™ Plates and the 3M™ Molecular Detection System, and is effective at neutralizing sanitizers commonly found in the food industry.

AOAC® Performance Tested MethodsSM Certificate 022104.*

Request a free trial at 3M.com/SampleCollection.

*The 3M Environmental Scrub Sampler with 10mL Wide Spectrum Neutralizer is the first and only sample collection device and neutralizing solution to be granted the AOAC® Performance Tested MethodSM certification.





Introducing a new class of capability.

You'll want to get to know the 3M™ Environmental Scrub Sampler with 10 mL Wide Spectrum Neutralizer, an innovative solution for environmental microbial sampling applications. It effectively neutralizes commonly used sanitizers and allows for up to a 96-hour sample hold time. The stick is designed to be ergonomic, metal detectable and features a convenient 1-step detachment for ease of use.

3M™ Environmental Scrub Sampler

- 1 Ergonomic stick design conforms to multiple grip types
- 2 Stick contains an additive for metal detectability
- 3 Bright color for visibility in food processing environments
- 4 1-step twist detachment
- 5 Maintains integrity when sampling rough surfaces
- 6 Multilayered non-woven composite product designed for pick-up, release and absorbency
- 7 Hydrated with a proprietary 10 mL Wide Spectrum Neutralizer
- 8 Scrub dot technology to disrupt biofilm and enhance sample collection



Description	Catalog #	Quantity
3M™ Environmental Scrub Sampler Stick with 10mL Wide Spectrum Neutralizer	ESS10WSN	100/case
3M™ Environmental Scrub Sampler with 10mL Wide Spectrum Neutralizer and Gloves	HES10WSN2G	100/case

Wide Spectrum Neutralizer

- Formulated to effectively neutralize common sanitizers, including quats, high acid cleaners, chlorine-based, hydrogen peroxide and peroxyacetic acid-based sanitizers
- Free from known allergenic components, known PCR-inhibitors, animal-derived materials (ADMs) and components made from genetically modified organisms (GMOs)
- Designed for compatibility with 3M™ Petrifilm™ Plates and the 3M™ Molecular Detection System
- Allows up to a 96-hour hold time

Distributed by:

NELSON JAMESON
INC.

800-826-8302 nelsonjameson.com