

3M Food Safety

3M™ Sample Collection Products



# Experience confidence

The complete 3M sample handling system

## Collect.

Better testing begins with better samples — and better samples are obtained when you use the best collection tools. 3M™ Sample Collection solutions provide quality products designed specifically for microbial surface sampling. Make environmental testing and product sampling easy and convenient with the full line of 3M Sample Collection products.

- Easy to use for product sampling and environmental testing
- Exceptional convenience and quality
- Increased productivity and efficiency
- Reliable and accurate — time after time
- Improved consistency from technician to technician and plant to plant



# Collect.

Sponges and swabs designed specifically for microbial surface sampling.

## 3M™ Sponge-Sticks\*

- Allows sampling without the need to directly handle the sponge, increasing operator efficiency with less risk of contamination
- Easier to reach ceilings, into drains and pipes, and around equipment
- Sturdy sponge snaps free from plastic handle after sampling
- Biocide-free cellulose sponge maintains *Listeria* organism viability within 72 hours of refrigeration\*\*
- Available in dry or pre-hydrated formats, with or without irradiated gloves



## 3M™ Hydrated Sponges

- Ready-to-use, pre-hydrated sponges for easy and convenient sample collection
- Available in a variety of pre-hydrated dilutions
- Biocide-free cellulose sponge maintains *Listeria* organism viability within 72 hours of refrigeration\*\*
- Durable sponge withstands scrubbing
- Sponge bag with write-on label
- Available with or without irradiated gloves



## 3M™ Dry-Sponges

- Dry sponge in a sample bag
- Ideal for those who prefer to use their own sampling broths or a broth not available with the 3M™ Hydrated Sponges
- Durable sponge made from biocide-free cellulose to help maintain organism viability



## 3M™ Extendable Poles

- Designed for use with the 3M™ Sponge Stick for easier sampling of remote and hard-to-reach areas
- Lightweight, sturdy aluminum construction
- Minimizes bending and stretching of plant personnel
- Safely sample ceilings and tanks



\*Formerly known as Spongesicle.

\*\*Technical data summary can be provided upon request.

| Catalog No.              | Description   | Qty. |
|--------------------------|---|------|
| <b>3M™ Sponge-Sticks</b> |   |      |
| SSL100                   | 3M™ Sponge-Stick, 1.5" x 3" Sponge, dry, 30 oz. bag   | 100  |
| SSL10BPW                 | 3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Buffered Peptone Water Broth, 30 oz. bag, 5 bags of 20  | 100  |
| SSL10LET                 | 3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Letheen Broth, 5 bags of 20                             | 100  |
| SSL10LET2G               | 3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Letheen Broth, 30 oz. bag, 2 gloves, 5 bags of 20       | 100  |
| SSL10DE                  | 3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL D/E Neutralizing Broth, 30 oz. bag, 5 bags of 20        | 100  |
| SSL10NB                  | 3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 30 oz. bag, 5 bags of 20           | 100  |
| SSL10NB2G                | 3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 30 oz. bag, 2 gloves, 5 bags of 20 | 100  |

|                             |   |     |
|-----------------------------|---|-----|
| <b>3M™ Hydrated Sponges</b> |   |     |
| HS10BPW                     | 3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Buffered Peptone Water Broth, 5 bags of 20           | 100 |
| HS10BPW2G                   | 3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Buffered Peptone Water Broth, 2 gloves, 5 bags of 20 | 100 |
| HS10DE2G                    | 3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL D/E Neutralizing Buffer, 2 gloves                    | 100 |
| HS10LET                     | 3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Letheen Broth, 5 bags of 20                          | 100 |
| HS10LET2G                   | 3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Letheen Broth, 2 gloves, 5 bags of 20                | 100 |
| HS119510NB                  | 3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 50 oz. bag                      | 100 |
| HS10NB                      | 3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 5 bags of 20                    | 100 |
| HS10NB2G                    | 3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 2 gloves, 5 bags of 20          | 100 |
| HS2410NB2G                  | 3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 24 oz. bag, 2 gloves            | 100 |

|                        |  |     |
|------------------------|--|-----|
| <b>3M™ Dry-Sponges</b> |  |     |
| BP135DS                | 3M™ Dry-Sponge w/4" String, 1.5" x 3" Sponge, 18 oz. bag | 100 |
| BP133ES                | 3M™ Dry-Sponge, 18 oz. printed bag                       | 100 |
| BP237SPG               | 3M™ Dry-Sponge, Sponge w/glove, 18 oz. bag               | 100 |

|                             |   |   |
|-----------------------------|---|---|
| <b>3M™ Extendable Poles</b> |   |   |
| XSLSSL10NB                  | 3M™ Extendable Pole Kit, (1) XSL-008 and (1) SSL-10NB | 1 |
| XSL008                      | 3M™ Extendable Pole w/Sponge Holder, 8 ft.            | 1 |



### 3M™ Quick Swabs

- Ready-to-use environmental swab quickly delivers 1mL of sample onto a 3M™ Petrifilm™ Plate
- Increased accuracy
- Improved consistency from technician to technician and plant to plant
- Labor savings in dilution mixing and pipetting
- Better productivity and efficiency
- Can be used wet or dry



### 3M™ Swab-Samplers

- Ideal for sampling hard-to-reach areas with long, pre-hydrated 9cm swab
- Premeasured fill volumes allow for qualitative or quantitative use
- Blue stick is easy to detect if dropped onto production area or floor
- Can stand alone or fit into a standard test tube rack
- Ample headspace for effective vortexing

### 3M™ Enviro Swabs

- Designed to minimize the risk of sample contamination and improve the detection of bacteria on all types of surfaces
- Large swab head with shaped tip to sample crevices for efficient collection
- Flexible plastic stem allows swabbing of difficult-to-access areas
- Free of substances that inhibit bacteria
- Specially formatted transport medium
- Easily re-sealed for transport



| Catalog No.            | Description    | Qty. |
|------------------------|----------------|------|
| <b>3M™ Quick Swabs</b> |                |      |
| 6432                   | 3M™ Quick Swab | 50   |
| 6433                   | 3M™ Quick Swab | 250  |

|                          |  |     |
|--------------------------|--|-----|
| <b>3M™ Swab-Samplers</b> |  |     |
| RS96010BPW               | 3M™ Swab-Sampler w/10mL Buffered Peptone Water Broth | 100 |
| RS96010DE                | 3M™ Swab-Sampler w/10mL D/E Neutralizing Broth       | 100 |
| RS9601LET                | 3M™ Swab-Sampler w/1mL Letheen Broth                 | 100 |
| RS9604LET                | 3M™ Swab-Sampler w/4mL Letheen Broth                 | 100 |
| RS9605LET                | 3M™ Swab-Sampler w/5mL Letheen Broth                 | 100 |
| RS96010LET               | 3M™ Swab-Sampler w/10mL Letheen Broth                | 100 |
| RS9604NB                 | 3M™ Swab-Sampler w/4mL Neutralizing Buffer           | 100 |
| RS96010NB                | 3M™ Swab-Sampler w/10mL Neutralizing Buffer          | 100 |

|                         |                      |    |
|-------------------------|----------------------|----|
| <b>3M™ Enviro Swabs</b> |                      |    |
| ENVSWB25                | 3M™ Enviro Swab      | 25 |
| ENVSWB25D               | 3M™ Enviro Swab, Dry | 25 |





