



Homogenizer Bags

Where quality and reliability are in the bag!



EASY TO OPEN AND EASY TO USE

- Available in both standard and filtered styles and sized to fit any 400 mL machine
- Write-on panel for labeling samples or placement of bar code label
- Filtered style allows for separation of the liquid extract from sample solids thereby preventing pipette blockage
- Filter with solid contents can be easily removed if desired
- Packaged in slide trays for easy access

STRONG ENOUGH TO ENDURE THE SAMPLE PREPARATION PROCESS, EVEN AT LOW TEMPERATURES


- Thick 3 mil polyethylene for toughness
- Additional additive for increased elasticity
- Extra wide 4 mm heat seal is designed to prevent leaking
- Filter section of the filtered bag is made of durable, perforated polyethylene film; the longer filter length gives extra reinforcement and support during the blending process

NO WORRIES ABOUT PRIOR CONTAMINATION OF THE BAG

- Gamma irradiated

BACKED BY 3M'S QUALITY REPUTATION

- Manufactured with proven experience and uncompromising quality



Ordering Information

Product	Catalog No.	Contents	Dimensions
3M Homogenizer Bag Standard 400	6468	500 bags/case	180x300mm
3M Homogenizer Bag Filtered 400	6469	200 bags/case	180x300mm

Additional Information

- Questions? U.S., call **1-800-328-6553**.
- To order, call **1-800-328-1671**.
- 3M Microbiology offers a full line of products to accomplish a variety of your microbial testing needs. For more product information, visit us at **www.3M.com/microbiology**.



3M Microbiology

3M Center, Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-800-228-3957
www.3M.com/microbiology
microbiology@mmm.com



Recycled Paper
40% pre-consumer
10% post-consumer

Printed in U.S.A.

© 3M 2003
70-2009-4616-1