

Thermo Scientific Automation Tips made to fit

Caliper/Zymark Liquid Handling Workstations

- Low percentage coefficient of variation
- Available in nonsterile, sterile or sterile filtered
- Validated on the actual workstation for which the tip is created

Workstations:

- Caliper Presto™
- Caliper Allegro™
- Caliper Zephyr™
- Caliper SciClone™
- Caliper/Zymark RapidPlate™



VOLUME: 50 µl

Cat. No.	Description	Qty
920-021	50 µl clear, sterile, filtered	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
920-261	50 µl clear, sterile	
920-262	50 µl clear, nonsterile	



VOLUME: 100 µl

Cat. No.	Description	Qty
921-021	80 µl clear, sterile, filtered	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
921-261	100 µl clear, sterile	
921-262	100 µl clear, nonsterile	



VOLUME: 200 µl

Cat. No.	Description	Qty
923-021	150 µl clear, sterile, filtered	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
923-261	200 µl clear, sterile	
923-262	200 µl clear, nonsterile	

Distributed by:

NELSON JAMESON
INC.

800-826-8302 nelsonjameson.com