

Accu-Chek Aviva Plus Test Strips



Features

- Accu-Chek Aviva Plus test strips, Easy-fill tip, the easy-fill tip quickly and easily absorbs blood sample
- Tiny sample Fills quickly with only a 0.6 µl drop of blood
- Wide sample area, The wide sample area simplifies blood application
- The Accu-Chek Aviva test strips can be used with the following Accu-Chek meters: Accu-Chek Aviva ,Accu-Chek Aviva Combo, Accu-Chek Aviva Connect , Accu-Chek Aviva Expert, Accu-Chek Aviva Insight

ACCU-CHEK® Aviva Plus

For testing glucose in fresh capillary whole blood
For single patient use only
For in vitro diagnostic use
Store between 36–86 °F

Use between 10–80 % relative humidity

Mutant variant of quinoprotein glucose dehydrogenase (Mut. Q-GDH, EC 1.1.5.2 from Acinetobacter spec.)

• **Fast 5-second test**

• **Finger, palm, or forearm testing***

Manufactured in the U.S.A. using U.S. and imported materials.

*Check with your healthcare professional before testing with blood from sites other than your fingertip. Do not use alternate site testing with the ACCU-CHEK Aviva Combo or the ACCU-CHEK Aviva Expert systems. Refer to the AST information in the package insert.

Cat. No. 06908268001
NDC No. 65702-0408-10

ACCU-CHEK® Aviva Plus

For questions, contact the ACCU-CHEK Customer Care Service Center toll-free at 1-800-858-8072. We offer assistance in many languages.

Si tiene preguntas, llame gratuitamente al ACCU-CHEK Customer Care Service Center al 1-800-858-8072. Ofrecemos asistencia en muchos idiomas.

ACCU-CHEK, ACCU-CHEK AVIVA, ACCU-CHEK AVIVA COMBO, and ACCU-CHEK AVIVA EXPERT are trademarks of Roche.

Roche Diagnostics
9115 Hague Road
Indianapolis, IN 46256

accu-chek.com

nelsonjameson.com

800-826-8302

NELSON JAMESON
INC.