

Product specification

Product:	PE plasters with detectable X ray pad
-----------------	---------------------------------------

1-Adhesive:

Adhesive type	Hotmelt adhesive
Coat weight (g/m2)	35-45
Adhesion power (N/25mm)	20±3
Loop Tack (N/25mm)	15±4
Shear (hrs)	≥20hrs

2- Wound pad:

Composition	gm/sqm	1- non adherent layer 100% PE 2- absorption layer 75% viscose 25% PP/PE 3-Backing : Aluminum sheet With a Barium sulphate ribbon 1 mm thickness and 4 mm width
TOTAL GRAMMAGE	g/m2	220+/-10%
Absorption capacity	%	≥400
Time of absorption	Sec	< 1

3- Release paper

Test	UNITS	TOLL.
GRAMMAGE	g/m2	75± 3



Distributed by:



800-826-8302 nelsonjameson.com

www.detectamet.co.uk

T: +44 (0)1759 304200 | F: +44 (0)1759 305236
E: sales@detectamet.com

Company Registration No. 5103699
VAT Registration No. GB 842 411357

Detectamet Ltd

Prospect House, 1 Halifax Way, Pocklington Ind. Estate
York, YO42 1NR, United Kingdom

4-Wrapping paper:

A) front side wrapping paper

TEST	UNITS	TOLL.
Total weight	g/m2	45+/-7

B) back side wrapping paper

Test	UNITS	TOLL.
Weight	g/m2	44+/-7

Sealing strength for the pouches = ≥ 1.0 N/15mm

Breathable and Hypoallergenic

Distributed by:

NELSON JAMESON INC.

800-826-8302 nelsonjameson.com

www.detectamet.co.uk

