



250ml Blood Bank Refrigerator



60ml



100ml



30ml

Block Heater style Thermometer



*Space Saving Dual Magnets
for quick release to record temperature.
Each magnet has an adhesive backing
to attach onto Stainless, ABS, or
Glass walls of refrigerators or freezers*

Environmentally Friendly Spirit Filled Thermometers (Non-Mercury)

CATALOG No.	Application	Range	Divisions	Height	Bottle Media
ACCBBS	Blood Banks Refrigerator	-5 to 20°C	0.5°C	145mm	Bio-Safe Glycol/Water - 250 ml
ACCB210S	Blood Banks Refrigerator	-2 to 10°C	0.1°C	180mm	Bio-Safe Glycol/Water - 250 ml
ACCF100S	Freezers	-25 to -5°C	0.1°C	120mm	Bio-Safe Glycol/Water - 60 ml
ACCF0101S	Freezers	-30 to 0°C	0.5°C	125mm	Bio-Safe Glycol/Water - 30 ml
ACCR20001S	Refrigerators	-2 to 10°C	0.1°C	240mm	Bio-Safe Glycol/Water - 60 ml
ACCR0201S	Refrigerators	-5 to 15°C	0.5°C	125mm	Bio-Safe Glycol/Water - 30 ml
ACCI300S	Incubators	15 to 30°C	0.1°C	180mm	Bio-Safe Glycol/Water - 60 ml
ACCI0301S	Incubators	18 to 50°C	0.5°C	135mm	Bio-Safe Glycol/Water - 30 ml
ACCI310S	Incubators	25 to 45°C	0.1°C	190mm	Bio-Safe Glycol/Water - 60 ml
ACCAM0401S	Ambient/Room	+10 to 30°C	0.5°C	140mm	Bio-Safe Glycol/Water - 30 ml
ACCOV070S	Oven	20 to 130°C	1.0°C	135mm	Sand - 30 ml
ACCOV700S	Oven	50 to 115°C	0.5°C	240mm	Sand - 60 ml
ACCOV080S	Oven	35 to 200°C	1.0°C	145mm	Sand - 100 ml Glass Bottle
ACCULF010S	Ultra-Low Freezers	-90 to 20°C	1.0°C	145mm	Sand - 60 ml
ACCCDB301	Block Heaters	18 to 60°C	0.5°C	125mm	Mineral Oil - 1 ml
ACCCDB40	Block Heaters	50 to 110°C	0.5°C	165mm	Mineral Oil - 1 ml

Mercury Filled Thermometers (Hg)

CATALOG No.	Application	Range	Divisions	Height	Bottle Media
ACCB210	Blood Banks Refrigerator	-2 to 10°C	0.1°C	180mm	Bio-Safe Glycol/Water - 250 ml
ACCF100	Freezers	-25 to -5°C	0.1°C	210mm	Bio-Safe Glycol/Water - 60 ml
ACCF0101	Freezers	-30 to 0°C	0.5°C	115mm	Bio-Safe Glycol/Water - 30 ml
ACCR200	Refrigerators	-2 to 10°C	0.1°C	180mm	Bio-Safe Glycol/Water - 60 ml
ACCR0201	Refrigerators	-5 to 15°C	0.5°C	125mm	Bio-Safe Glycol/Water - 30 ml
ACCI300	Incubators	15 to 30°C	0.1°C	180mm	Bio-Safe Glycol/Water - 60 ml
ACCI0301	Incubators	18 to 50°C	0.5°C	135mm	Bio-Safe Glycol/Water - 30 ml
ACCI310	Incubators	25 to 45°C	0.1°C	190mm	Bio-Safe Glycol/Water - 60 ml
ACCOV070	Oven	20 to 130°C	1.0°C	135mm	Sand 30ml
ACCOV700	Oven	95 to 115°C	0.1°C	200mm	Sand 60ml
ACCOV080	Oven	35 to 200°C	1.0°C	145mm	Sand 100ml