

# Aquachek®






Water Quality  
Test Strips for

**High-Range  
Chlorine**

100 Test Strips  
Cat. 28902-00

2013LB 8/04



**DIRECTIONS:** Dip a test strip into your water sample and remove immediately. Hold strip level, pad side up, for 30 seconds. Compare the pad to the Free Chlorine color chart on the label.


**IMPORTANT: KEEP CAP ON TIGHT BETWEEN USES. STORE AT ROOM TEMPERATURE.**

**DIRECCIONES:** Sumerja una cinta de análisis en su muestra del agua y retírela inmediatamente. Sostenga la cinta nivelada, con las almohadillas hacia arriba, durante 30 segundos. Compare la almohadilla de la cinta con la carta de colores de cloro libre en la etiqueta.

**IMPORTANTE: MANTENGA LA TAPA BIEN APRETADA ENTRE UN USO Y OTRO. ALMACENE EL PRODUCTO A TEMPERATURA AMBIENTE.**

**DIRECTIONS :** Trempez une bandelette dans votre échantillon d'eau et retirez-la immédiatement. Tenez la bandelette à l'horizontale (carrés vers le haut) pendant 30 secondes. Comparez le carré de la bandelette au tableau des couleurs du chlore libre, sur l'étiquette.

**IMPORTANT : BIEN REFERMER APRÈS UTILISATION. CONSERVER A TEMPERATURE AMBIANTE.**



USE BY DATE  
ON BOTTOM

Hach Company, P.O. Box 389, Loveland, CO 80539 U.S.A.  
(800) 227-4224 Outside U.S.A. (970) 669-3050