




The Right Choice for Bacteria Detection



Because your environmental testing results can be significantly influenced by the swab you use, it is important to choose your swab carefully. The 3M™ Enviro Swab is designed to minimize risk of sample contamination and improve the detection of bacteria on all types of surfaces.

- Large swab head with shaped tip and flexible stem
- Free of substances inhibiting bacteria
- Specially formulated transport medium
- Easily re-sealed for transport

3M™ Enviro Swab

► Improved Detection Rate

- Large swab head with shaped tip to sample crevices for efficient collection
- Flexible plastic stem allows swabbing of difficult-to-access areas
- Foam head is free of substances that inhibit bacteria
- Specially formulated transport medium maximizes survival of damaged or stressed microorganisms and neutralizes effect of sanitizing agents

► Convenient, Cost-Effective and Saves Time



- No need to purchase and prepare transport medium
- Swab density and strength suitable for numerous applications
- Pre-enrichment and incubation steps can be performed in the tube
- Easy to re-seal for transport
- Dry swabs also available

Order Information

► 3M™ Enviro Swab

3M Product	Description	Catalog Number	Quantity
3M™ Enviro Swab	Pre-enriched	ENVSWB25	25/pack
	Dry	ENVSWB25D	25/pack

In trials, the 3M™ Enviro Swab gave superior recovery of *Listeria* and *Salmonella* when compared with ordinary cotton swabs.



3M Food Safety offers a full line of products and solutions that work together for consistent, reliable results. For more information, contact your local 3M representative or visit us at www.3M.com/foodsafety.



3M Food Safety

3M Center
Building 275-5W-05
St. Paul, MN 55144-1000
USA
1-800-328-6553
www.3M.com/foodsafety

3M is a trademark of 3M. Please recycle.
Printed in U.S.A. © 3M 2011. All rights reserved. 70-2009-9501-0