



Transport-Swab | Sample Collection

Transport Swab is the ideal sampling device for customers looking for an environmental sample collection device containing more than 1mL of pre-filled diluent. The device is available in 4mL & 10mL volumes. The prefilled bottles come sterilized and with ample headspace for additional sample to be added. The Rayon swab tip sits in the broth, allowing organisms to release from the swab bud. Transport Swab makes sample collection efficient and easy. Collect your samples by swabbing the surface, recapping, and the liquid is ready to be transferred to the media of our choice.




Find support documents, instructional videos, and more at www.hygiena.com

Benefits:

- Ready to use and convenient
- All-in-one sample collection system
- More broth for high-volume samples
- Perfect for plates or dried media film
- Available broths: Letheen & Buffered Peptone Water
- 18 month shelf life at refrigerated temperatures (2-8°C)

Transport Swab Procedure:



Catalog No.	Description	Quantity
HTS-4	Letheen 4 ml	100
HTS-10	Letheen 10 ml	100
HTS-10BPW	Buffered Peptone Water 10 ml	100

Distributed by:

NELSON JAMESON INC.

800-826-8302 nelsonjameson.com

REV A 112013