

# 3M™ Petrifilm™ *Salmonella* Confirmation Disks

## Product Information

MEDIUM:	3M Petrifilm <i>Salmonella</i> Confirmation Disks (6538/6539)				
ISO/BSI CERTIFICATE OF REGISTRATION NUMBER:	Manufactured at Brookings, SD, USA (ISO 9001:2008, FM 14552)				
DATE OF EXPIRATION/ LOT NUMBER:	Expiry and lot number indicated on each package. Lot number indicated on each packaged disk.				
FORMULATION:	Differential substrate				
METHOD OF PREPARATION:	Differential substrate coated onto films. For use, place into 3M Petrifilm <i>Salmonella</i> Express Plates. See product instructions for detailed instructions.				
EFFICACY CHECK:	<p>The disks are tested for the appropriate functioning of the differential substrate. The following test determines the function:</p> <table><tr><td><u>Enzyme</u></td><td><u>Acceptable Result</u></td></tr><tr><td>Proprietary Differential Enzyme</td><td>Blue color development where enzyme is placed within five minutes</td></tr></table>	<u>Enzyme</u>	<u>Acceptable Result</u>	Proprietary Differential Enzyme	Blue color development where enzyme is placed within five minutes
<u>Enzyme</u>	<u>Acceptable Result</u>				
Proprietary Differential Enzyme	Blue color development where enzyme is placed within five minutes				
PACKAGING:	<p>Pack size: 5 individually packaged disks per foil pouch Film grade: Plastic foil laminate Seal integrity check: Delamination seal integrity test performed</p>				
MEASUREMENT SYSTEMS CALIBRATION AND TRACEABILITY:	<p>Incubator temperature 3M internal calibration Minimum calibration once per year for all equipment</p>				
MEDIA QUALITY STATEMENT:	Certificate of Analysis available on <a href="http://3m.com/foodsafety">3m.com/foodsafety</a>				
SHELF LIFE:	30 months from date of manufacture				
STORAGE CONDITIONS:	Store at temperatures 2 to 8°C				

Distributed by:



800-826-8302 [nelsonjameson.com](http://nelsonjameson.com)