


PRODUCT INFORMATION

Item # 485-3352

 **ATAGO U.S.A., Inc.**



MASTER-53T

MASTER-53T / MASTER-53PT analog refractometer are capable of measuring juice, coffee, beverages and also processed food such as sauce, ketchup, low calorie jam, etc.

The MASTER-T / PT series is a good alternative choice when the water resistant feature is not a priority.

The MASTER-53T is made of metal, and the [MASTER-53PT](#) is plastic.

For salty and acidic samples, choose the plastic body, [MASTER-53PT](#).

Specifications

| | |
|---------------|---|
| Model | MASTER-53T |
| Cat.No. | 2352 |
| Range | Brix : 0.0 to 53.0% |
| Minimum Scale | 0.2% |
| Accuracy | ±0.2% |
| Repeatability | ±0.1% |
| Optional | <ul style="list-style-type: none">• SUPER EXTRACTOR : RE-29401• 10% Sucrose Solution (±0.03%) : RE-110010• 20% Sucrose Solution (±0.03%) : RE-110020• 30% Sucrose Solution (±0.03%) : RE-110030• 40% Sucrose Solution (±0.04%) : RE-110040• 50% Sucrose Solution (±0.05%) : RE-110050• Daylight Plate for Small Volume Samplings : RE-2311-87M• Daylight Plate : RE-2315-60M• MASTER Carrying Case: RE-29410* Calibration Certificate: NIST-traceable certificate. This will increase the lead time of the order by 2 to 3 weeks. Contact an ATAGO representative for details. |

NELSON JAMESON
INC.

nelsonjameson.com

800-826-8302