

Nasco SAMPLING

901 Janesville Avenue, Fort Atkinson, WI 53538-0901

Nasco Carcass Sampling Kit Material Specification Sheet B01324WA

The carcass sampling kit is made of the following components:

1. B01245WA-Whirl-Pak® Speci-Sponge® Bag

Size: 4 ½" W x 9" L (11 ½ x 23 cm).

Capacity: 18 oz. (532 ml).

Film: Blend of low-linear and low-density virgin polyethylene, 2.5 mil (.064 mm) thick, FDA approved for food contact.

Write-on Strip: Horizontal band of thermoset white ink printed on film, compliant with CONEG resolution for heavy metals.

Closing Tab: Aluminum wires covered with plastic tape by a patented process that also covers the wire ends, reduces bag puncture, forms leak-proof seal when closed.

Pull Tabs: Made of Tyvek®.

Sponge: contains 1 ½" x 3" (3.8 x 7.6 cm) dehydrated, sterile, natural cellulose sponge, tested to be non-inhibitory.

Sterilization: Sterilized with ethylene oxide gas.

ISO Certification: Manufactured under a quality system certified to **ISO9001 & ISO13485**.

Packaging: 25 bags/kit.

2. Swabbing Templates & Gloves

--Two paper templates with cut out area measuring 10 x 10 cm for sampling beef and pork, and 5 x 10 cm for sampling turkeys.

--One pair of polyethylene gloves.

--Templates and gloves are packaged together in a sealed bag and sterilized, 25 packs/kit.

3. Buffered Peptone Water

--25 ml of buffered peptone water in a PET (polyethylene terephthalate) bottle, gamma irradiated, 25 bottles/kit.

--PH @ 25C-7.2, +/-0.2

--Shelf life is 12 months from date of manufacture. Expiration date is marked on the outside of each box.

WEB NascoSampling.com

PHONE 1.800.558.9595

FAX 1.800.372.1236

EMAIL contact@nascosampling.com

Distributed by:

NELSON JAMESON
INC.

800-826-8302 nelsonjameson.com