


Code 2934

QAC Dual Range  **LaMotte**

50 Test Strips

Insta-Test[®] ANALYTIC



Using a cup-size sample, immerse strip for **2 seconds. SHAKE ONCE TO REMOVE EXCESS WATER.** Wait **15 seconds** and **immediately compare** to color chart.



- Keep wet fingers out of the bottle.
- Cap bottle tightly.
- Store in a cool, dry place.



LaMotte Company
Chestertown, MD 21620
<http://www.lamotte.com>

Lot# 3429 Exp. 9/15

