


AlerTox® Sticks



Rapid tests for allergen detection in
food, drinks, and surfaces for internal
Quality Control.



Crustacean



Egg



Hazelnut



Fish



β-Lactoglobulin



Almond



Mustard Seed



Casein



Peanut



Soy



Total Milk


AlerTox® Sticks advantages:

- » Quick and Reliable
- » No Cross-Reactivity
- » Results in 10 minutes
- » High Sensitivity

AlerTox® Sticks are a practical tool for the control of allergen content in raw materials, final products, working surfaces.

Industrial application

- » Semi-quantification of allergens in food
- » Food quality control
- » Trace allergen contamination in food, liquid and surfaces
- » Implementation and fulfillment of HACCP, ISO 22000, BRC, and IFS compliance



Product Name	Ordering No.	Reference No.	Quantity
AlerTox® Sticks Crustacean	KIT3036	KT-6126	10 tests
AlerTox® Sticks Egg	KIT3025	KT-5899	10 tests
AlerTox® Sticks Hazelnut	KIT3035	KT-6129	10 tests
AlerTox® Sticks Fish	KIT3038	KT-6131	10 tests
AlerTox® Sticks β-Lactoglobulin	KIT3018	KT-5782	10 tests
AlerTox® Sticks Almond	KIT3033	KT-6127	10 tests
AlerTox® Sticks Mustard Seed	KIT3040	KT-6370	10 tests
AlerTox® Sticks Casein	KIT3021	KT-5781	10 tests
AlerTox® Sticks Peanut	KIT3031	KT-6128	10 tests
AlerTox® Sticks Soy	KIT3028	KT-6125	10 tests
AlerTox® Sticks Total Milk	KIT3015	KT-6495	10 tests

NOTE: As with any qualitative or semiquantitative test, results should be periodically validated with a quantitative test performed in an outside laboratory.

