






Page & Pedersen  
International, Ltd.

# Comparison: QuickCheck™ Cryoscopes



## SPECIFICATIONS

	<b>QuickCheck™ Flex</b>		<b>QuickCheck™ Tech</b>		<b>QuickCheck™ IV/GR QuickCheck™ Select</b>	
						
Lab or Field Use	Portable w/Carrying Cover		Lab		Lab	
Calibration Channels	1	Yes	1	Yes	1	3
Raw Milk	Yes		Yes		Yes	Yes
Lactose Hydrolysis	No		No		Yes	Yes
Ice Cream	No		No		No	Yes
Freezing Point Range	0 to -1000m <sup>0</sup> C or m <sup>0</sup> H		0 to -1000m <sup>0</sup> C or m <sup>0</sup> H		0 to -1000m <sup>0</sup> C or m <sup>0</sup> H	
					Channel A Milk: -600m <sup>0</sup> C Channel B Lactose Hydrolysis: -600 to -992m <sup>0</sup> C Channel C Ice Cream Mix -1 <sup>0</sup> C to -2.5 <sup>0</sup> C	
DC Power Input	Yes 11-14V		No		No	
AC Power Input	W/Optional External Power Supply		Yes 100 - 260V, 50-60Hz		Yes 100 - 260V, 50-60Hz	
Test Time	~2 min.		90 sec - 2 min		90 sec - 2 min	
Read Out Units m <sup>0</sup> c or m <sup>0</sup> H	Yes		Yes		Yes	
Ambient Air Operating Temperature	15-40 <sup>0</sup> C 59-104 <sup>0</sup> F		15-38 <sup>0</sup> C 59-100 <sup>0</sup> F		15-38 <sup>0</sup> C 59-100 <sup>0</sup> F	
Printer	External (Optional)		External (Optional)		External (Optional)	
Connector	PC/USB to memory Stick		PC/USB to memory Stick		PC/USB to memory Stick	
Graphic Display	No		No		Yes	
Digital Readout	Yes		Yes		Yes	

DEDICATED ANALYTICAL

158 West Main Street, Hopkinton, MA 01748 U.S.A.

TOOLS FOR THE DAIRY SUPPLY CHAIN

© (508) 435 5966, (508) 435 8198

info@pagepedersen.com

www.pagepedersen.com