



Bartovation LLC.
2483 47th St.
Astoria, NY 11103

Product Technical Fact Sheet

Part Number: PQT03V50

Description: High level Quick Response QAC Test Strips 0-1500 PPM

Application:

The 0-1500 ppm QAC test strips are used for testing disinfection/sanitizing solutions of quaternary ammonium mixtures in water. The test strip and color chart have been calibrated for use with Oasis 146 (Ecolab). Oasis 146 concentrate consists of the following active ingredients:

| | |
|---|-------|
| Alkyl dimethyl benzyl ammonium chloride | 3.00% |
| Octyl decyl dimethyl ammonium chloride | 2.25% |
| Dioctyl dimethyl ammonium chloride | 0.90% |
| Didecyl dimethyl ammonium chloride | 1.35% |

The test strip can be used with other similar Quat formulas, however, the specific application should be verified before use.

Test Range: 0-1500 ppm

Test Increments: Color chart calibrated at 0, 200, 400, 750, 1000 and 1500 ppm using Oasis 146

Accuracy: +/- ½ color chart unit

Detection Limit: 50 ppm

Storage Recommendations: Store in original packaging in a cool (20-30C), dry, place out of direct sunlight.

Shelf-Life: Two years from date of manufacture when stored properly in original packaging.

Interferences:

Formula additives can adversely influence the results. Use with formulations other than Oasis 146 should be verified first.

Temperature can adversely affect the results. Temperatures of 90-110F can cause strips to read up to ½ unit high. Temperatures above 110F can shift the results 1 unit high. Testing at room temperature is recommended.

Hard water (above 500 ppm) can cause the strips to read high, up to 1 unit.

Instructions for Use:

1. Dip test strip into the solution to be tested for 1-2 seconds.
2. Remove test strip and immediately compare to color chart.

Chemistry Behind the Test:

The QAC test strips use various surfactants, buffers, and indicators to employ the protein-error of indicators effect.

Distributed by:
NELSON JAMESON
INC.
800-826-8302 nelsonjameson.com