



Bartovation LLC.
2483 47th St.
Astoria NY 11103
contact@bartovation.com

Product Technical Fact Sheet

Part Number: PCL08V100

Description: Chlorine Test Paper, 10-200 ppm [Vial of 100 Paper Strips]

Application:

The 0-200 ppm Chlorine test paper is used for testing sanitization and disinfecting levels of total available chlorine in water. The Chlorine test papers can be used in facilities where food preparation areas, equipment, and general disinfection are strictly controlled.

Test Range: 0-200 ppm

Test Increments: Color chart calibrated at 0, 10, 50, 100, and 200 ppm using sodium hypochlorite

Accuracy: +/- 1/2 color chart unit

Detection Limit: 5 ppm

Storage Recommendations: Store in original packaging in a cool (20-30C), dry, place out of direct sunlight.

Shelf-Life: Two years from date of manufacture when stored properly in original packaging.

Interferences: Test strips may react to other oxidizing materials, peracetic acid, peroxide, etc.

Instructions for Use:

1. Use dry fingers to remove a chlorine test paper strip from the vial.
2. Dip the test strip into the solution to be tested, and do not agitate.
3. Remove and compare to the color chart immediately.

Chemistry Behind the Test:

The test paper is based on the Potassium Iodide Starch reaction. The intensity of the color developed is proportional to the approximate amount of available (total) chlorine present.

Other Information:

The Chlorine test papers were developed to work with solutions prepared from concentrated solutions of sodium hypochlorite.

Disclaimer:

The test papers are not intended for use to diagnose, treat, or monitor any medical condition.

Distributed by:



800-826-8302 nelsonjameson.com